





Infant Incubator

Model: 111



Infant Incubator 111 Specification

- Digital Air Temperature Thermometer**
 - Resolution 0.1°C
 - Controllable range 20 to 37°C Override to 39°C
- Electrical requirements**
 - Power Supply 220V 50Hz 3A
- Alarms**
 - Air Sensor fault
 - Skin Sensor fault
 - Air circulation failure
 - Power failure
 - Aire set deviation +/- 2° C
 - Low Battery
 - Low Water
- Digital Skin Temperature thermometer**
 - Resolution 0.1°C
 - Controllable range 28 to 37.5°C
- Digital Humidity Sensor**
 - Resolution 1 % RH
 - Controllable range 30 - 90 % RH
- Heater - Power indicator**
 - LED Power
- Performance**
 - Internal noise level less than 48 dB
 - Warm-up time 30 minutes
 - Power consumption 500 W
 - Humidification up to 90% RH
 - Disposable filter
- Dimension and weight**
 - Mattress 34 * 60 Cm
 - Length, width, height 92, 63, 145
 - Weight 90 Kg
- Available Option**
 - Iv pole
 - Drawer
 - Height adjustment