




Manufacture & Supply of Neonatal System

**توسان Tosan**

# Intensive Phototherapy

Model 222,224,225



Manufacture & Supply of Neonatal System

**توسان Tosan**

## Intensive Phototherapy

### Performance

- LED Tube 8, 12, 16 (9W)
- Code of LEDs T8-600B02
- Touch screen Control panel
- Digital Display
- Digital's display to measure inside temperature
- Digital's display to measure Baby temperature
- Time and Date display
- Timer for measuring the effective life of the lamps
- Treatment time measuring
- Adjustable LED lights 50% or 100%
- Exposure angle 360 °
- Inside ventilation fans for ventilation of internal air and cooling the electrical parts

### Alarms


- Fan failure Audible and visible
- End of treatment time Visible
- Power failure Visible
- Overheat Audible and visible

### Electrical Requirements

- Power Supply 110 – 240 V AC, 50 Hz- 1.5 A

### Dimensions & weight

- Lenght 71 Cm
- Width 74 Cm
- Height 121 Cm
- Weight 37 Kg



Head Office

First south floor, Chehelsotoon Bldg.  
North Allameh St. Saadatabad, Tehran-Iran

Tel.: +98 21 22097416-7

Fax: +98 21 22362814

E-mail: [info@tosanco.com](mailto:info@tosanco.com)

Web: [www.tosanco.com](http://www.tosanco.com)